
Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | PKC alpha Antibody |
| Catalog Number: | R12-3330 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, ELISA |
| Dilutions: | PKC alpha Antibody |
| Modification: | Unmodified/Total |
| Source: | Rabbit Polyclonal |
| Purification: | The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human PKC alpha. |
| Concentration: | 1 mg/ml |
| Gene Name: | PRKCA |
| Other Name: | kinase PKC-alpha; KPCA; PKC III; PKC-A; PKC-alpha; PKC-III; PRKACA; PRKCA; Protein kinase C alpha type |
| Human Swiss Prot No.: | P17252 |
| Formulation: | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Specificity: | PKC alpha Antibody antibody detects endogenous levels of PKC alpha. |
| Storage and Stability: | Stable at -20°C for at least 1 year. |