

Antibodies Product Datasheet

Product Name:	RPR1A Antibody
Catalog Number:	C51194
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C-terminal RPR1A
Concentration:	1 mg/ml
Gene Name:	RPRD1A P15RS
Protein Name:	Regulation of nuclear pre-mRNA domain-containing protein 1A (Cyclin-dependent kinase inhibitor 2B-related protein) (p15INK4B-related protein)
Other Name:	Regulation of nuclear pre-mRNA domain-containing protein 1A (Cyclin-dependent kinase inhibitor 2B-related protein) (p15INK4B-related protein)
Human Swiss Prot No.:	Q96P16
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RPR1A at Human, Mouse
Storage and Stability:	Store at -20° C/1 year