

Antibodies Product Datasheet

Product Name:	PSMC6 Antibody
Catalog Number:	C14017
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PSMC6. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	PSMC6
Protein Name:	26S protease regulatory subunit 10B
Other Name:	PSMC6; SUG2; 26S protease regulatory subunit 10B; 26S proteasome AAA-ATPase subunit RPT4; Proteasome 26S subunit ATPase 6; Proteasome subunit p42
Human Swiss Prot No.:	P62333
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PSMC6 Polyclonal Antibody detects endogenous levels of PSMC6 protein.
Storage and Stability:	-20°C/1 year