

## Antibodies Product Datasheet

<b>Product Name:</b>	eIF2B $\delta$ Antibody
<b>Catalog Number:</b>	C19274
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human EIF2B4. AA range:226-275
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	EIF2B4
<b>Protein Name:</b>	Translation initiation factor eIF-2B subunit delta
<b>Other Name:</b>	EIF2B4; EIF2BD; Translation initiation factor eIF-2B subunit delta; eIF-2B GDP-GTP exchange factor subunit delta
<b>Human Swiss Prot No.:</b>	Q9UI10
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	eIF2B $\delta$ Polyclonal Antibody detects endogenous levels of eIF2B $\delta$ protein.
<b>Storage and Stability:</b>	-20°C/1 year