

Antibodies Product Datasheet

Product Name:	TADA2L Antibody
Catalog Number:	C10181
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ADA2L. AA range:206-255
Concentration:	1 mg/ml
Gene Name:	TADA2A
Protein Name:	Transcriptional adapter 2-alpha
Other Name:	TADA2A; TADA2L; KL04P; Transcriptional adapter 2-alpha; Transcriptional adapter 2-like; ADA2-like protein
Human Swiss Prot No.:	O75478
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TADA2L Polyclonal Antibody detects endogenous levels of TADA2L protein.
Storage and Stability:	-20°C/1 year