

Antibodies Product Datasheet

Product Name:	CNOT2 Antibody
Catalog Number:	B8387
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CNOT2. AA range:67-116
Concentration:	1 mg/ml
Gene Name:	CNOT2
Protein Name:	CCR4-NOT transcription complex subunit 2
Other Name:	CNOT2; CDC36; NOT2; HSPC131; MSTP046; CCR4-NOT transcription complex subunit 2; CCR4-associated factor 2
Human Swiss Prot No.:	Q9NZN8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CNOT2 Polyclonal Antibody detects endogenous levels of CNOT2 protein.
Storage and Stability:	-20°C/1 year