

Antibodies Product Datasheet

Product Name:	TSC22D1 Antibody
Catalog Number:	R12-3653
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	TSC22D1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TSC22D1.
Concentration:	1 mg/ml
Gene Name:	TSC22D1
Other Name:	TSC22 domain family protein 1; Transforming growth factor beta-1-induced transcript 4 protein; Regulatory protein TSC-22; TGFB-stimulated clone 22 homolog; Cerebral protein 2; TSC22D1; KIAA1994; TGFB114; TSC22
Human Swiss Prot No.:	Q15714
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	TSC22D1 Antibodyantibody detects endogenous levels of TSC22D1.
Storage and Stability:	Stable at -20°C for at least 1 year.