

Antibodies Product Datasheet

Product Name:	CATB (Cleaved-Leu80) Antibody
Catalog Number:	L0458
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CATB (Cleaved-Leu80)
Concentration:	1 mg/ml
Gene Name:	CTSB CPSB
Protein Name:	CATB (Cleaved-Leu80)
Other Name:	Cathepsin B (EC 3.4.22.1) (APP secretase) (APPS) (Cathepsin B1) [Cleaved into: Cathepsin B light chain; Cathepsin B heavy chain]
Human Swiss Prot No.:	P07858
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human CATB (Cleaved-Leu80, protein was cleaved amino acid sequence between 79-80)
Storage and Stability:	-20°C/1 year