

Antibodies Product Datasheet

Product Name:	Adducin β Antibody
Catalog Number:	C12041
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ADD2. AA range:471-520
Concentration:	1 mg/ml
Gene Name:	ADD2
Protein Name:	Beta-adducin
Other Name:	ADD2; ADDB; Beta-adducin; Erythrocyte adducin subunit beta
Human Swiss Prot No.:	P35612
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Adducin β Polyclonal Antibody detects endogenous levels of Adducin β protein.
Storage and Stability:	-20°C/1 year