

Antibodies Product Datasheet

Product Name:	CXXC5 Antibody
Catalog Number:	C48295
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CXXC5
Concentration:	1 mg/ml
Gene Name:	CXXC5 HSPC195 TCCCIA00297
Protein Name:	CXXC-type zinc finger protein 5 (CF5) (Putative MAPK-activating protein PM08) (Putative NF-kappa-B-activating protein 102) (Retinoid-inducible nuclear factor) (RINF)
Other Name:	
Human Swiss Prot No.:	Q7LFL8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CXXC5 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year