

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | CaMKIV Antibody |
| Catalog Number: | B0831 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, IHC-p, IF/ICC, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human CaMK4. AA range:166-215 |
| Concentration: | 1 mg/ml |
| Gene Name: | CAMK4 |
| Protein Name: | Calcium/calmodulin-dependent protein kinase type IV |
| Other Name: | CAMK4; CAMK; CAMK-GR; CAMKIV; Calcium/calmodulin-dependent protein kinase type IV; CaMK IV; CaM kinase-GR |
| Human Swiss Prot No.: | Q16566 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | CaMKIV Polyclonal Antibody detects endogenous levels of CaMKIV protein. |
| Storage and Stability: | -20°C/1 year |