

Antibodies Product Datasheet

Product Name:	Hox-A1/B1/D1 Antibody
Catalog Number:	C10372
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HOXA1/B1/D1. AA range:216-265
Concentration:	1 mg/ml
Gene Name:	HOXB1
Protein Name:	Homeobox protein Hox-B1
Other Name:	HOXB1; HOX2I; Homeobox protein Hox-B1; Homeobox protein Hox-2I; HOXA1; HOX1F; Homeobox protein Hox-A1; Homeobox protein Hox-1F; HOXD1; HOX4; HOX4G; Homeobox protein Hox-D1; Homeobox protein Hox-GG
Human Swiss Prot No.:	P14653/P49639/Q9GZZ0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Hox-A1/B1/D1 Polyclonal Antibody detects endogenous levels of Hox-A1/B1/D1 protein.
Storage and Stability:	-20°C/1 year