

Antibodies Product Datasheet

Product Name:	Amyloid- β Antibody
Catalog Number:	C30857
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Amyloid- β at AA range: 631-680
Concentration:	1 mg/ml
Gene Name:	APP
Protein Name:	Amyloid beta A4 protein, Amyloid- β , A β
Other Name:	Amyloid beta A4 protein (ABPP) (APPI) (APP) (Alzheimer disease amyloid protein) (Cerebral vascular amyloid peptide) (CVAP) (PreA4) (Protease nexin-II) (PN-II) [Cleaved into: N-APP; Soluble APP-alpha (S-APP-alpha); Soluble APP-beta (S-APP-beta); C99; Beta-P05067
Human Swiss Prot No.:	P05067
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Amyloid- β Polyclonal Antibody detects endogenous levels of Amyloid- β
Storage and Stability:	-20°C/1 year