

Antibodies Product Datasheet

Product Name:	CSTF2T Antibody
Catalog Number:	R12-2658
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	CSTF2T Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CSTF2T.
Concentration:	1 mg/ml
Gene Name:	CSTF2T
Other Name:	Cleavage stimulation factor subunit 2 tau variant; CF-1 64 kDa subunit tau variant; Cleavage stimulation factor 64 kDa subunit tau variant; CSTF 64 kDa subunit tau variant; TauCstF-64; CSTF2T; KIAA0689
Human Swiss Prot No.:	Q9H0L4
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	CSTF2T Antibodyantibody detects endogenous levels of CSTF2T.
Storage and Stability:	Stable at -20°C for at least 1 year.