

Antibodies Product Datasheet

Product Name:	LDOC1L Antibody
Catalog Number:	C11219
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LDOC1L. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	LDOC1L
Protein Name:	Protein LDOC1L
Other Name:	LDOC1L; MAR6; MART6; Protein LDOC1L; Leucine zipper protein down-regulated in cancer cells-like; Mammalian retrotransposon-derived protein 6
Human Swiss Prot No.:	Q6ICC9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LDOC1L Polyclonal Antibody detects endogenous levels of LDOC1L protein.
Storage and Stability:	-20°C/1 year