

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | MMP-23 Antibody |
| Catalog Number: | C0271 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human MMP-23. AA range:341-390 |
| Concentration: | 1 mg/ml |
| Gene Name: | MMP23A |
| Protein Name: | Matrix metalloproteinase-23 |
| Other Name: | MMP23A; MMP21; MMP23B; MMP21; MMP22; Matrix metalloproteinase-23; MMP-23; Femalysin; MIFR-1; Matrix metalloproteinase-21; MMP-21; Matrix metalloproteinase-22; MMP-22 |
| Human Swiss Prot No.: | O75900 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | MMP-23 Polyclonal Antibody detects endogenous levels of MMP-23 protein. |
| Storage and Stability: | -20°C/1 year |