

Antibodies Product Datasheet

Product Name:	ALK (Phospho-Tyr1586) Antibody
Catalog Number:	A8456
Reactivity:	Human:Y1586, Mouse:Y1592
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ALK (Phospho-Tyr1586)
Concentration:	1 mg/ml
Gene Name:	ALK
Protein Name:	ALK (Phospho-Tyr1586)
Other Name:	ALK tyrosine kinase receptor (EC 2.7.10.1) (Anaplastic lymphoma kinase) (CD antigen CD246)
Human Swiss Prot No.:	Q9UM73
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of ALK (Phospho-Tyr1586) at Human:Y1586, Mouse:Y1592
Storage and Stability:	-20°C/1 year