

Antibodies Product Datasheet

Product Name:	JIP3 Antibody
Catalog Number:	R12-2941
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	JIP3 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human JIP3.
Concentration:	1 mg/ml
Gene Name:	MAPK8IP3
Other Name:	C-jun-amino-terminal kinase interacting protein 3; JIP-3; JIP3; JNK MAP kinase scaffold protein 3; JNK-interacting protein 3; JNK/SAPK-associated protein 1; JSAP1; KIAA1066; MAPK8IP3; Mitogen- activated protein kinase 8-interacting protein 3; SYD2; Sunday Q9UPT6
Human Swiss Prot No.:	Q9UPT6
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	JIP3 Antibodyantibody detects endogenous levels of JIP3.
Storage and Stability:	Stable at -20°C for at least 1 year.