

Antibodies Product Datasheet

Product Name:	GNL3L Antibody
Catalog Number:	R12-2819
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	GNL3L Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GNL3L.
Concentration:	1 mg/ml
Gene Name:	GNL3L
Other Name:	Guanine nucleotide-binding protein-like 3-like protein; GNL3L
Human Swiss Prot No.:	Q9NVN8
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	GNL3L Antibodyantibody detects endogenous levels of GNL3L.
Storage and Stability:	Stable at -20°C for at least 1 year.