
Antibodies Product Datasheet

Product Name:	Pax-5 Antibody
Catalog Number:	C30750
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human PAX5. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	PAX5
Protein Name:	Paired box protein Pax-5
Other Name:	PAX5; Paired box protein Pax-5; B-cell-specific transcription factor; BSAP
Human Swiss Prot No.:	Q02548
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Pax-5 Polyclonal Antibody detects endogenous levels of Pax-5 protein.
Storage and Stability:	-20°C/1 year