

## Antibodies Product Datasheet

<b>Product Name:</b>	XIAP Antibody
<b>Catalog Number:</b>	B0599
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human XIAP. AA range:53-102
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	XIAP
<b>Protein Name:</b>	E3 ubiquitin-protein ligase XIAP
<b>Other Name:</b>	XIAP; API3; BIRC4; IAP3; E3 ubiquitin-protein ligase XIAP; Baculoviral IAP repeat-containing protein 4; IAP-like protein; ILP; hILP; Inhibitor of apoptosis protein 3; IAP-3; hIAP-3; hIAP3; X-linked inhibitor of apoptosis protein; X-linked I
<b>Human Swiss Prot No.:</b>	P98170
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	XIAP Polyclonal Antibody detects endogenous levels of XIAP protein.
<b>Storage and Stability:</b>	-20°C/1 year