

Antibodies Product Datasheet

Product Name:	WNT2B Antibody
Catalog Number:	C51190
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human N-terminal WNT2B
Concentration:	1 mg/ml
Gene Name:	WNT2B WNT13
Protein Name:	Protein Wnt-2b (Protein Wnt-13)
Other Name:	Protein Wnt-2b (Protein Wnt-13)
Human Swiss Prot No.:	Q93097
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of WNT2B at Human, Mouse
Storage and Stability:	Store at -20° C/1 year