

Antibodies Product Datasheet

Product Name:	GDF-15 Antibody
Catalog Number:	C15982
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GDF15. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	GDF15
Protein Name:	Growth/differentiation factor 15
Other Name:	GDF15; MIC1; PDF; PLAB; PTGFB; Growth/differentiation factor 15; GDF-15; Macrophage inhibitory cytokine 1; MIC-1; NSAID-activated gene 1 protein; NAG-1; NSAID-regulated gene 1 protein; NRG-1; Placental TGF-beta; Placental bone morphogenetic
Human Swiss Prot No.:	Q99988
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GDF-15 Polyclonal Antibody detects endogenous levels of GDF-15 protein.
Storage and Stability:	-20°C/1 year