

Antibodies Product Datasheet

Product Name:	RDH11 Antibody
Catalog Number:	C49813
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human RDH11
Concentration:	1 mg/ml
Gene Name:	RDH11 ARSDR1 PSDR1 CGI-82
Protein Name:	Retinol dehydrogenase 11 (EC 1.1.1.300) (Androgen-regulated short-chain dehydrogenase/reductase 1) (HCV core-binding protein HCBP12) (Prostate short-chain dehydrogenase/reductase 1) (Retinal reductase 1) (RalR1)
Other Name:	
Human Swiss Prot No.:	Q8TC12
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RDH11 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year