

Antibodies Product Datasheet

Product Name:	UBE2E1 Antibody
Catalog Number:	C50902
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human UBE2E1
Concentration:	1 mg/ml
Gene Name:	UBE2E1 UBCH6
Protein Name:	Ubiquitin-conjugating enzyme E2 E1 (EC 6.3.2.19) (UbcH6) (Ubiquitin carrier protein E1) (Ubiquitin-protein ligase E1)
Other Name:	
Human Swiss Prot No.:	P51965
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of UBE2E1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year