
Antibodies Product Datasheet

Product Name:	Ribosomal Protein L27A Antibody
Catalog Number:	C14165
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RPL27A. AA range:99-148
Concentration:	1 mg/ml
Gene Name:	RPL27A
Protein Name:	60S ribosomal protein L27a
Other Name:	RPL27A; 60S ribosomal protein L27a
Human Swiss Prot No.:	P46776
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ribosomal Protein L27A Polyclonal Antibody detects endogenous levels of Ribosomal Protein L27A protein.
Storage and Stability:	-20°C/1 year