

Antibodies Product Datasheet

Product Name:	S10A8 Antibody
Catalog Number:	C46052
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	S100A8 CAGA CFAG MRP8
Protein Name:	Protein S100-A8 (Calgranulin-A) (Calprotectin L1L subunit) (Cystic fibrosis antigen) (CFAG) (Leukocyte L1 complex light chain) (Migration inhibitory factor-related protein 8) (MRP-8) (p8) (S100 calcium-binding protein A8) (Urinary stone protein band A)
Other Name:	
Human Swiss Prot No.:	P05109
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	S10A8 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year