

Antibodies Product Datasheet

Product Name:	P44/42 MAPK (ERK1/2) Monoclonal Antibody(6B1)
Catalog Number:	N1553
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:50-100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of P44/42 MAPK (ERK1/2)
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	
Other Name:	MAPK3; ERK1; PRKM3; Mitogen-activated protein kinase 3; MAP kinase 3; MAPK 3; ERT2; Extracellular signal-regulated kinase 1; ERK-1; Insulin-stimulated MAP2 kinase; MAP kinase isoform p44; p44-MAPK; Microtubule-associated protein 2 kinase; p
Human Swiss Prot No.:	P27361/P28482
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	P44/42 MAPK(ERK1/2) protein detects endogenous levels of P44/42 MAPK(ERK1/2)
Storage and Stability:	-20°C/1 year