

Antibodies Product Datasheet

Product Name:	PERP Antibody
Catalog Number:	C51082
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PERP
Concentration:	1 mg/ml
Gene Name:	PERP KCP1 KRTCAP1 PIGPC1 THW
Protein Name:	p53 apoptosis effector related to PMP-22 (Keratinocyte-associated protein 1) (KCP-1) (P53-induced protein PIGPC1) (Transmembrane protein THW)
Other Name:	
Human Swiss Prot No.:	Q96FX8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PERP at Human, Mouse
Storage and Stability:	Store at -20° C/1 year