

Antibodies Product Datasheet

Product Name:	NHE-8 Antibody
Catalog Number:	C18840
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p:1:50-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC9A8. AA range:451-500
Concentration:	1 mg/ml
Gene Name:	SLC9A8
Protein Name:	Sodium/hydrogen exchanger 8
Other Name:	SLC9A8; KIAA0939; NHE8; Sodium/hydrogen exchanger 8; Na(+)/H(+) exchanger 8; NHE-8; Solute carrier family 9 member 8
Human Swiss Prot No.:	Q9Y2E8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NHE-8 Polyclonal Antibody detects endogenous levels of NHE-8 protein.
Storage and Stability:	-20°C/1 year