

Antibodies Product Datasheet

Product Name:	DRP-2 Antibody
Catalog Number:	R12-2714
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	DRP-2 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DRP-2.
Concentration:	1 mg/ml
Gene Name:	DPYSL2
Other Name:	Collapsin response mediator protein 2; CRMP2; Dihydropyrimidinase related protein-2; DPYL2; DPYSL2; DRP-2; N2A3; TOAD-64; Turned on after division, 64 kDa protein; ULIP 2 protein; ULIP2
Human Swiss Prot No.:	Q16555
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	DRP-2 Antibody antibody detects endogenous levels of DRP-2.
Storage and Stability:	Stable at -20°C for at least 1 year.