

Antibodies Product Datasheet

Product Name:	RPAP2 Antibody
Catalog Number:	C18453
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RPAP2. AA range:460-509
Concentration:	1 mg/ml
Gene Name:	RPAP2
Protein Name:	Putative RNA polymerase II subunit B1 CTD phosphatase RPAP2
Other Name:	RPAP2; C1orf82; Putative RNA polymerase II subunit B1 CTD phosphatase RPAP2; RNA polymerase II-associated protein 2
Human Swiss Prot No.:	Q8IXW5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RPAP2 Polyclonal Antibody detects endogenous levels of RPAP2 protein.
Storage and Stability:	-20°C/1 year