

Antibodies Product Datasheet

Product Name:	Hrs (Phospho-Tyr216) Antibody
Catalog Number:	P12-1027
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human Hrs (phospho Tyr216)
Concentration:	1 mg/ml
Gene Name:	HGS
Protein Name:	Hepatocyte growth factor-regulated tyrosine kinase substrate
Other Name:	HGS; HRS; Hepatocyte growth factor-regulated tyrosine kinase substrate; Hrs; Protein pp110
Human Swiss Prot No.:	O14964
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Hrs (Y216) Polyclonal Antibody detects endogenous levels of Hrs protein only when phosphorylated at Y216.
Storage and Stability:	-20°C/1 year