
Antibodies Product Datasheet

Product Name:	Mast cell tryptase mouse mAb(PT2025)
Catalog Number:	N1748
Reactivity:	Human
Applications:	IHC, WB
Dilutions:	IHC, WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Mast cell tryptase
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	Mast cell tryptase
Other Name:	
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Mast cell tryptase
Storage and Stability:	-20°C/1 year