

Antibodies Product Datasheet

Product Name:	GW182 Antibody
Catalog Number:	C51066
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GW182
Concentration:	1 mg/ml
Gene Name:	TNRC6A CAGH26 KIAA1460 TNRC6
Protein Name:	Trinucleotide repeat-containing gene 6A protein (CAG repeat protein 26) (EMSY interactor protein) (GW182 autoantigen) (Protein GW1) (Glycine-tryptophan protein of 182 kDa)
Other Name:	
Human Swiss Prot No.:	Q8NDV7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GW182 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year