

## Antibodies Product Datasheet

<b>Product Name:</b>	Glucosidase II $\beta$ Antibody
<b>Catalog Number:</b>	C10365
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human GLU2B. AA range:81-130
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PRKCSH
<b>Protein Name:</b>	Glucosidase 2 subunit beta
<b>Other Name:</b>	PRKCSH; G19P1; Glucosidase 2 subunit beta; 80K-H protein; Glucosidase II subunit beta; Protein kinase C substrate 60.1 kDa protein heavy chain; PKCSH
<b>Human Swiss Prot No.:</b>	P14314
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Glucosidase II $\beta$ Polyclonal Antibody detects endogenous levels of Glucosidase II $\beta$ protein.
<b>Storage and Stability:</b>	-20°C/1 year