

Antibodies Product Datasheet

Product Name:	MEPCE Antibody
Catalog Number:	C48332
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MEPCE
Concentration:	1 mg/ml
Gene Name:	MEPCE BCDIN3
Protein Name:	7SK snRNA methylphosphate capping enzyme (MePCE) (EC 2.1.1.-) (Bicoid-interacting protein 3 homolog) (Bin3 homolog)
Other Name:	
Human Swiss Prot No.:	Q7L2J0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MEPCE at Human, Mouse
Storage and Stability:	Store at -20° C/1 year