

Antibodies Product Datasheet

Product Name:	MAGE-1 Antibody
Catalog Number:	C0254
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAGE-1. AA range:260-309
Concentration:	1 mg/ml
Gene Name:	MAGEA1
Protein Name:	Melanoma-associated antigen 1
Other Name:	MAGEA1; MAGE1; MAGE1A; Melanoma-associated antigen 1; Antigen MZ2-E; Cancer/testis antigen 1.1; CT1.1; MAGE-1 antigen
Human Swiss Prot No.:	P43355
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MAGE-1 Polyclonal Antibody detects endogenous levels of MAGE-1 protein.
Storage and Stability:	-20°C/1 year