

Antibodies Product Datasheet

Product Name:	Anti-Prostatic Acid Phosphatase (PSAP) [ABT-PSAP]
Catalog Number:	V0118
Clone ID:	ABT-PSAP
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	P15309
Gene Name:	ACPP
Other Name:	Prostatic acid phosphatase;PAP;5'-NT;Ecto-5'-nucleotidase;TMPase;Thiamine monophosphatase
Isotype	Mouse IgG
Localization:	N/A
Positive Tissue Control:	N/A
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	Anti-Prostatic Acid Phosphatase (PSAP) [ABT-PSAP] detects endogenous levels of total Prostatic Acid Phosphatase (PSAP) protein.
Storage and Stability:	-20°C/1 year