

Antibodies Product Datasheet

Product Name:	KDM6A/UTX Antibody
Catalog Number:	C51055
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KDM6A/UTX
Concentration:	1 mg/ml
Gene Name:	KDM6A UTX
Protein Name:	Lysine-specific demethylase 6A (EC 1.14.11.-) (Histone demethylase UTX) (Ubiquitously-transcribed TPR protein on the X chromosome) (Ubiquitously-transcribed X chromosome tetratricopeptide repeat protein)
Other Name:	
Human Swiss Prot No.:	O15550
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KDM6A/UTX at Human, Mouse
Storage and Stability:	Store at -20° C/1 year