
Antibodies Product Datasheet

Product Name:	Desmin Monoclonal Antibody(1B12)
Catalog Number:	N1550
Reactivity:	Human,Rat,Mouse
Applications:	IF/ICC,IHC-p
Dilutions:	IF: 1:50-200 IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Desmin
Concentration:	1 mg/ml
Gene Name:	DES
Protein Name:	Desmin
Other Name:	DES; Desmin
Human Swiss Prot No.:	P17661
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Desmin protein detects endogenous levels of Desmin
Storage and Stability:	-20°C/1 year