

Antibodies Product Datasheet

Product Name:	L-Selectin Antibody
Catalog Number:	C22131
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from L-Selectin . at AA range: 60-140
Concentration:	1 mg/ml
Gene Name:	SELL
Protein Name:	L-selectin
Other Name:	SELL; LNHR; LYAM1; L-selectin; CD62 antigen-like family member L; Leukocyte adhesion molecule 1; LAM-1; Leukocyte surface antigen Leu-8; Leukocyte-endothelial cell adhesion molecule 1; LECAM1; Lymph node homing receptor; TQ1; gp90-MEL; CD62L
Human Swiss Prot No.:	P14151
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	L-Selectin Polyclonal Antibody detects endogenous levels of L-Selectin protein.
Storage and Stability:	-20°C/1 year