

## Antibodies Product Datasheet

<b>Product Name:</b>	ALDH3A1 Antibody
<b>Catalog Number:</b>	R12-2017
<b>Reactivity:</b>	Human, Rat
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human ALDH3A1. AA range:236-285
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ALDH3A1
<b>Protein Name:</b>	Aldehyde dehydrogenase dimeric NADP-preferring
<b>Other Name:</b>	ALDH3A1; ALDH3; Aldehyde dehydrogenase, dimeric NADP-preferring; ALDHIII; Aldehyde dehydrogenase 3; Aldehyde dehydrogenase family 3 member A1
<b>Human Swiss Prot No.:</b>	P30838
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	ALDH3A1 Polyclonal Antibody detects endogenous levels of ALDH3A1 protein.
<b>Storage and Stability:</b>	-20°C/1 year