

## Antibodies Product Datasheet

<b>Product Name:</b>	c-Abl (Phospho-Thr735) Antibody
<b>Catalog Number:</b>	A9438
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:1000-2000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized phospho peptide around human c-Abl (Thr735)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ABL1 ABL JTK7
<b>Protein Name:</b>	c-Abl (Thr735)
<b>Other Name:</b>	Tyrosine-protein kinase ABL1 (EC 2.7.10.2) (Abelson murine leukemia viral oncogene homolog 1) (Abelson tyrosine-protein kinase 1) (Proto-oncogene c-Abl) (p150)
<b>Human Swiss Prot No.:</b>	P00519
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Human c-Abl (phospho-Thr735)
<b>Storage and Stability:</b>	-20°C/1 year