

Antibodies Product Datasheet

Product Name:	ICK (Phospho-Tyr159) Antibody
Catalog Number:	P12-1118
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	ICK (Phospho-Tyr159) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human ICK around the phosphorylation site of Tyr159.
Concentration:	1 mg/ml
Gene Name:	ICK
Other Name:	Intestinal cell kinase; hICK; MAK-related kinase; MRK; Laryngeal cancer kinase 2; LCK2
Human Swiss Prot No.:	Q9UPZ9
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ICK (Phospho-Tyr159) Antibody detects endogenous levels of ICK only when phosphorylated at Tyr159.
Storage and Stability:	Stable at -20°C for at least 1 year.