

Antibodies Product Datasheet

Product Name:	Secretory phospholipase A2 Antibody
Catalog Number:	C51068
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Secretory phospholipase A2
Concentration:	1 mg/ml
Gene Name:	PLA2G2A PLA2B PLA2L RASF-A
Protein Name:	Phospholipase A2, membrane associated (EC 3.1.1.4) (GIIC sPLA2) (Group IIA phospholipase A2) (Non-pancreatic secretory phospholipase A2) (NPS-PLA2) (Phosphatidylcholine 2-acylhydrolase 2A)
Other Name:	
Human Swiss Prot No.:	P14555
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Secretory phospholipase A2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year