

Antibodies Product Datasheet

Product Name:	PDGFR- β Monoclonal Antibody
Catalog Number:	N0511
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human PDGFR- β expressed in E. Coli.
Concentration:	
Gene Name:	PDGFRB
Protein Name:	Beta-type platelet-derived growth factor receptor
Other Name:	PDGFRB; PDGFR; PDGFR1; Platelet-derived growth factor receptor beta; PDGF-R-beta; PDGFR-beta; Beta platelet-derived growth factor receptor; Beta-type platelet-derived growth factor receptor; CD140 antigen-like family member B; Platelet-deri
Human Swiss Prot No.:	P09619
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	PDGFR- β Monoclonal Antibody detects endogenous levels of PDGFR- β protein.
Storage and Stability:	-20°C/1 year