

Antibodies Product Datasheet

Product Name:	NOTCH3 (Cleaved-Val1662) Antibody
Catalog Number:	L0444
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NOTCH3 (Cleaved-Val1662)
Concentration:	1 mg/ml
Gene Name:	NOTCH3
Protein Name:	NOTCH3 (Cleaved-Val1662)
Other Name:	Neurogenic locus notch homolog protein 3 (Notch 3) [Cleaved into: Notch 3 extracellular truncation; Notch 3 intracellular domain]
Human Swiss Prot No.:	Q9UM47
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat NOTCH3 (Cleaved-Val1662, protein was cleaved amino acid sequence between 1661-1662)
Storage and Stability:	-20°C/1 year