

Antibodies Product Datasheet

Product Name:	FGF-6 Antibody
Catalog Number:	C30078
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human FGF6. AA range:159-208
Concentration:	1 mg/ml
Gene Name:	FGF6
Protein Name:	Fibroblast growth factor 6
Other Name:	FGF6; HST2; HSTF2; Fibroblast growth factor 6; FGF-6; Heparin secretory-transforming protein 2; HST-2; HSTF-2; Heparin-binding growth factor 6; HBGF-6
Human Swiss Prot No.:	P10767
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FGF-6 Polyclonal Antibody detects endogenous levels of FGF-6 protein.
Storage and Stability:	-20°C/1 year