

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Chk1 Monoclonal Antibody
<b>Catalog Number:</b>	N0152
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant fragment of Chk1 expressed in E. Coli.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CHK1
<b>Protein Name:</b>	Serine/threonine-protein kinase Chk1
<b>Other Name:</b>	CHEK1; CHK1; Serine/threonine-protein kinase Chk1; CHK1 checkpoint homolog; Cell cycle checkpoint kinase; Checkpoint kinase-1
<b>Human Swiss Prot No.:</b>	O14757
<b>Formulation:</b>	Purified antibody in PBS containing 0.03% sodium azide.
<b>Specificity:</b>	Chk1 Monoclonal Antibody detects endogenous levels of Chk1 protein.
<b>Storage and Stability:</b>	-20°C/1 year